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## TRANSFORMATION / DISSOLUTION TESTING – TEST SYSTEM AND SPECIATION ANALYSIS

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
Figure 1: An example showing that the dissolution of a test substance (metal compound) is pH-dependent.

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
Figure 2: Example of the concentration of a dissolved metal measured over time (test performed in triplicate vessels).

F4

F5

## Background and aims

The European REACH and CLP regulations will inevitably require the submission of data concerning the ecotoxicological potential of numerous metals and metal compounds. Therefore, it may be necessary to test such compounds according to Chemicals Testing Monograph No. 29: Guidance Document on Transformation / Dissolution of Metals and Metal Compounds in Aqueous Media (2001).

## Approach

The guidance document describes a test system that determines the rate at which and the extent to which metals and sparingly-soluble metal compounds can produce soluble, available ionic and other metal-bearing species in aqueous media. The test conditions should be representative of the aqueous environment, and the test media are therefore based on reconstituted water according to ISO 6341 with a pH range of 5.5 to 8.5. The guidance document also demands accurate and comprehensible procedures (e.g. exact loadings, temperature monitoring, and mass balance analysis at the test endpoint). A system was developed that fulfills these test requirements for the transformation / dissolution of various metal substances and also redox-sensitive metals. The tests were carried out in accordance with the requirements of good laboratory practice (GLP). We developed and implemented methods adapted for particular needs for the quantitation of dissolved metal ions and species by ICP-OES or ICP-MS. These methods were validated by measuring certified reference materials, fortified samples, reagent and method blanks and recalibration samples along with the test samples.

## Results

Successful tests were carried out on a number of metal compounds. As expected, dissolution tended to be clearly dependent on pH (Fig. 1) and saturation was often observed during testing (Fig. 2).

The available analytical instruments are advantageous because they allow multi-element analysis if required, high sample throughput, low limits of quantification and a wide linear dynamic range. Different (redox) species may be formed during the transformation / dissolution tests and these must be stabilized, and can then be separated by coupling liquid chromatography to ICP-MS (Fig. 3).


Figure 3: Example of liquid chromatography used to separate two metal species from the OECD 29 test, followed by detection using ICP-MS.

## Contact / Ansprechpartner

Dr. Heinz Rüdell

Tel: +49 2972 302 - 301

heinz.ruedel@ime.fraunhofer.de

Dr. Thorsten Klawonn

Tel: +49 2972 302 - 119

thorsten.klawonn@ime.fraunhofer.de

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Consultants and industry

Figure 4: Liquid chromatography coupled to ICP-MS.

Figure 5: Transformation/dissolution testing at Fraunhofer IME.